## Phytopathologia Mediterranea

The international journal of the Mediterranean Phytopathological Union



**Citation:** S. Palavouzis, A. Triantafyllopoulou, C. Tsoukas, A.K. Tzima, E. Paplomatas (2023) Dual labelled probe assays for differentiation of *Botryosphaeria dothidea, Neofusicoccum mediterraneum* and *Neofusicoccum parvum*, based on polymorphisms in the MAT1-2-1 gene. *Phytopathologia Mediterranea*62(2):177-182.doi:10.36253/ phyto-14375

Accepted: April 24, 2023

Published: July 23, 2023

**Copyright:** © 2023 S. Palavouzis, A. Triantafyllopoulou, C. Tsoukas, A.K. Tzima, E. Paplomatas. This is an open access, peer-reviewed article published by Firenze University Press (http://www.fupress.com/pm) and distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

**Data Availability Statement:** All relevant data are within the paper and its Supporting Information files.

**Competing Interests:** The Author(s) declare(s) no conflict of interest.

**Editor:** José R. Úrbez Torres, Agriculture and Agri-Food Canada, Summerland, British Columbia, Canada.

### ORCID:

SP: 0000-0002-2305-4293 AT: 0000-0002-9460-0466 CT: 0000-0002-9460-0466 AKT: 0000-0002-1178-8222 EP: 0000-0002-2929-0676 Short Notes - 12th Special issue on Grapevine Trunk Diseases

# Dual labelled probe assays for differentiation of *Botryosphaeria dothidea*, *Neofusicoccum mediterraneum* and *Neofusicoccum parvum*, based on polymorphisms in the MAT1-2-1 gene

STAVROS PALAVOUZIS, ALEXANDRA TRIANTAFYLLOPOULOU, CHRISTOS TSOUKAS, ALIKI K. TZIMA\*, EPAMINONDAS PAPLOMATAS\*

Laboratory of Phytopathology, Department of Crop Science, School of Plant Sciences, Agricultural University of Athens, Greece

\*Corresponding author. E-mail: epaplom@aua.gr; aliki@aua.gr

**Summary.** Botryosphaeriaceous fungi are widespread, and cause serious diseases in many economically important crops. *Botryosphaeria dothidea, Neofusicoccum mediterraneum* and *N. parvum* are the most important members of this family in the Mediterranean region. These fungi are frequently isolated from the same host, which together with their extensive and increasing host range necessitates development of rapid and reliable diagnostic tools. Species boundaries within the *Botryosphaeriaceae* have been defined based on phylogenetic analyses of multiple gene sequences, including those of mating type genes. The MAT1-2-1 gene displayed high sequence variability between Botryosphaeriaceous species, so was selected as the target for development of a definitive diagnostic tool. This paper outlines a new and robust molecular tool, composed of three TaqMan assays based on polymorphisms located in the MAT1-2-1 gene of *B. dothidea, N. mediterraneum* and *N. parvum*. Each assay differentiated the target species from other *Botryosphaeriaceae*, and from non-target fungi.

Keywords. Diagnosis, mating type genes, molecular marker, real time PCR.

Fungal pathogens in the *Botryosphaeriaceae* have been reported as the causal agents of severe losses in a wide range of hosts, including agricultural, horticultural and forest plants (Slippers *et al.*, 2007; Vakalounakis *et al.*, 2019; Aiello *et al.*, 2020; Guarnaccia *et al.*, 2020; Batista *et al.*, 2021; Luna *et al.*, 2022; Guarnaccia *et al.*, 2023). Symptoms caused by these pathogens include cankers and dieback in twigs and branches, leaf necroses, blight, and fruit rots (Chen et al., 2014; Marsberg et al., 2017). Important hosts of these pathogens are pistachio (*Pistacia vera* L.), grapevine (*Vitis vinifera*) and citrus. The main species causing high yield losses in Greece and other Mediterranean countries, where pistachio and grapevine are cultivated, are *Botryosphaeria dothidea*, *Diplodia seriata*, *Neofusicoccum vitifusiforme*, *N. mediterraneum* and *N. parvum* (Lazzizera *et al.*, 2008; Chen *et al.*, 2015; Stempien *et al.*, 2017; Bezerra *et al.*, 2021). Members of the *Botryosphaeriaceae* are known to reproduce sexually, through ascospores formed in perithecia, and asexually from conidia produced in pycnidia. *Botryosphaeria dothidea* and *N. parvum* are homothallic, bearing the MAT1-1-1 and MAT1-2-1 genes in single strains, whereas *N. vitifusiforme*, and *N. mediterraneum* are heterothallic. For *N. vitifusiforme*, only MAT1-1-1 strains have been reported, whereas for *N. mediterraneum*, strains bearing MAT1-1-1 or MAT1-2-1 genes have been described (Lopes *et al.*, 2017; Marsberg *et al.*, 2017).

Species differentiation within the Botryosphaeriaceae has been mainly based on phylogenetic analyses of multiple gene sequences, such as the internal transcribed spacer region (ITS), the elongation factor 1-alpha, (TEF1-a), the beta tubulin, and the RNA polymerase II subunit (RPB2) (Pavlic et al., 2009; Chen et al., 2015). Although an RFLP analysis tool based on the ITS region could distinguish the main Botryosphaeriaceous species, this could not differentiate five species, including N. parvum (Slippers et al., 2007). More rapid molecular tools, such as PCR primers, have been developed for detection of Botryosphaeriaceous species only at the genus level, or for differentiation of Botryosphaeria dothidea from Neofusicoccum spp. (Ridgway et al., 2011; Palavouzis et al., 2022). Pathogen phylogenies based on mating type genes were comparable to multigene analyses for species differentiation. Notably, the MAT1-2-1 gene displayed high sequence variability between Botryosphaeriaceous species (Lopes et al., 2017).

With many species causing disease in an increasing range of hosts (Vakalounakis *et al.*, 2019; Batista *et al.*, 2020), rapid and accurate detection at the species level is important for the effective management of these pathogens, including quarantine measures and epidemiological studies. The present study reports development of robust and dual labelled probe molecular tools ("TaqMan" technology), targeting the MAT1-2-1 gene, for differentiation among *B. dothidea*, *N. mediterraneum* (MAT1-2-1 strains) and *N. parvum*, from fungal cultures and *in planta*.

Initially, gene sequences of MAT1-2-1 from N. mediterraneum, N. parvum and B. dothidea retrieved from the NCBI database, and sequences for each species derived from isolates of our own collection (deposited in GenBank under the Accession Numbers OQ632937 for B. dothidea, OQ632936 for N. parvum and OQ596433 for N. mediterraneum), were aligned using Clustal Omega (Madeira et al., 2019). Differences between sequences of the three species were detected using the MEGA software (Kumar et al., 2018). Based on sequence polymorphisms, dual labelled probes and primers were designed for MAT1-2-1 genes of B. dothidea (homothallic), N. mediterraneum MAT1-2-1 strain and N. parvum (homothallic). Primers were designed so that nucleotide polymorphisms specific for the target sequence were positioned at the 3' end of each primer, and PCR amplicons were smaller than 200 bp. Dual labelled probes were designed to have a Tm 8-10°C higher than the primers with a sequence length of 15-30 bp, avoiding a G nucleotide at the 5' end of the probe so that quenching of the 5' fluorophore was prevented. Nucleotide polymorphisms conferring probe specificity were placed close to the middle of the probe sequence. Primers and probes were screened for self-dimers, heterodimers and hairpins, using primer 3 plus software. To check target specificity, blast search against non-target Diplodia seriata MAT 1-2-1 sequences showed low homology with designed primers and probes. Primers and probes were then synthesized by eurofins genomics, labelling the 5' end of all probes (for B. dothidea, N. mediterraneum and N. parvum) with fluorescein (FAM) and the 3' end with the Black Hole Quencher-1 (BHQ1) (Table 1).

**Table 1.** List of primers and probe sequences for TaqMan qPCR assays, developed for detection and quantification of *Botryosphaeria dothidea, Neofusicoccum mediterraneum* and *Neofusicoccum parvum*, based on polymorphisms in the corresponding MAT1-2-1 genes. Polymorphisms of primers and probes between the different species are indicated with nucleotides letters in bold font.

Primer/probe name	Sequence (5'->3')	Target Gene	Product size (bp)	
Botdo2-1_295F	TCGCATCCTCTTCCCTCC <b>TG</b>			
Botdo2-1_380R	AGGCCAAGACCTGCTGA <b>AGT</b>	Botryosphaeria dothidea MAT1-2-1	86	
Tqm-Botdo2-1_316pb	ATACGTCGCACCCGCTCCCAAC	WIA11-2-1		
Neomed2-1_260F	GTCCGCGCTCCAGTCA <b>TC</b>	Neofusicoccum		
Neomed2-1_373R	AGGCTGAGGAGTGGAAACC	mediterraneum MAT1-	114	
Tqm-Neomed2-1_324pb	CGACCCTCCATGCTGACGGCG	2-1		
Neopar2-1_292F	CTGACCTTGTCCAGCACG	Neofusicoccum parvum		
Neopar2-1_424F	1_424F GCTGAGAAGCCGAA <b>GGTG</b>		133	
Tqm-Neopar2-1_377pb	CGAACTTCCGCGCCAATGGTATCAAC	MAT1-2-1		


MAT1-2-1 gene target	Template	Mean Ct	St. Dev. (3 reps)	
B. dothidea	B. dothidea 8 ng μL <sup>-1</sup>	20.6	0.2	
	<i>B. dothidea</i> 2 ng $\mu$ L <sup>-1</sup>	21.1	0.2	
	<i>B. dothidea</i> 200 pg μL <sup>-1</sup>	24.6		
	B. dothidea 200 pg $\mu L^{\text{-1}}$ & Pistacia vera 8 ng $\mu L^{\text{-1}}$	23.5	0.1	
	B. dothidea 20 pg $\mu L^{-1}$	28.1	0.3	
	<i>B. dothidea</i> 2 pg $\mu$ L <sup>-1</sup>	31.3	0.3	
	<i>B. dothidea</i> 200 fg $\mu$ L <sup>-1</sup>	35.7	0.3	
	<i>N. parvum</i> 8 ng $\mu$ L <sup>-1</sup>	35.6	1.1	
	N. mediterraneum "58" (MAT1-2-1) 8 ng $\mu$ L <sup>-1</sup>	undetermined		
	N. mediterraneum "28" (MAT1-1-1) 8 ng μL <sup>-1</sup>	undetermined		
	<i>P. chlamydospora</i> 8 ng μL <sup>-1</sup>	undetermined		
	Diplodia seriata 8 ng µL <sup>-1</sup>	undetermined		
	Pistacia vera 8 ng µL <sup>-1</sup>	undetermined / 35.5	NA	
N. mediterraneum	N. mediterraneum "58" (MAT1-2-1) 8 ng µL <sup>-1</sup>	21.7	0.1	
	N. mediterraneum "28" (MAT1-1-1) 8 ng $\mu$ L <sup>-1</sup>	undetermined		
	<i>N. mediterraneum</i> 2 ng $\mu$ L <sup>-1</sup>	24.7	0.0	
	N. mediterraneum 200 pg $\mu$ L <sup>-1</sup> 28.3	0.3		
	N. mediterraneum 200 pg $\mu L^{\text{-}1}$ & Pistacia vera 8 ng $\mu L^{\text{-}1}$	26.0	0.3	
	<i>N. mediterraneum</i> 20 pg $\mu$ L <sup>-1</sup>	31.9	0.3	
	<i>N. mediterraneum</i> 2 pg $\mu$ L <sup>-1</sup>	35.2	1.2	
	<i>B. dothidea</i> 8 ng $\mu$ L <sup>-1</sup>	36.5	0.4	
	<i>N. parvum</i> 8 ng $\mu$ L <sup>-1</sup>	35.8	0.7	
	<i>P. chlamydospora</i> 8 ng μL <sup>-1</sup>	undetermined		
	Diplodia seriata 8 n µL <sup>-1</sup>	undetermined		
	Pistacia vera 8 ng µL-1	37.5	1.2	
N. parvum	Diplodia seriata 8 n μL <sup>-1</sup> undetermined           Pistacia vera 8 ng μL <sup>-1</sup> 37.5           rvum         N. parvum 8 ng μL <sup>-1</sup> 22.4	0.3		
B. dothidea 8 ng $\mu$ L <sup>-1</sup> 36.5         N. parvum 8 ng $\mu$ L <sup>-1</sup> 35.8         P. chlamydospora 8 ng $\mu$ L <sup>-1</sup> undetermined         Diplodia seriata 8 n $\mu$ L <sup>-1</sup> undetermined         Pistacia vera 8 ng $\mu$ L <sup>-1</sup> 37.5         N. parvum       N. parvum 8 ng $\mu$ L <sup>-1</sup> 22.4         N. parvum 2 ng $\mu$ L <sup>-1</sup> 24.2         N. parvum 200 pg $\mu$ L <sup>-1</sup> & Pistacia vera 8 ng/ $\mu$ L <sup>-1</sup> 27.9         N. parvum 200 pg $\mu$ L <sup>-1</sup> & Pistacia vera 8 ng/ $\mu$ L <sup>-1</sup> 27.7	0.1			
	0.1			
	N. parvum 200 pg μL <sup>-1</sup> & Pistacia vera 8 ng/ μL <sup>-1</sup>	27.7	1.0	
	<i>N. parvum</i> 20 pg $\mu$ L <sup>-1</sup>	31.3	0.4	
	<i>N. parvum</i> 2 pg $\mu$ L <sup>-1</sup>	34.5	0.7	
	<i>B. dothidea</i> 8 ng $\mu$ L <sup>-1</sup>	36.3	0.5	
	N. mediterraneum "58" (MAT1-2-1) 8 ng $\mu$ L <sup>-1</sup>	undetermined		
	N. mediterraneum "28" (MAT1-1-1) 8 ng $\mu$ L <sup>-1</sup>	undetermined		
	<i>P. chlamydospora</i> 8 ng μL <sup>-1</sup>	undetermined		
	Diplodia seriata 8 n µL <sup>-1</sup>	undetermined		
	<i>Pistacia vera</i> 8 ng $\mu$ L <sup>-1</sup>	36.0	0.1	

**Table 2.** Mean Ct values (and standard deviations) for different fungi in three dual labelled probe assays targeting the MAT1-2-1 gene for diagnosis of *Botryosphaeria dothidea, Neofusicoccum mediterraneum* (MAT 1-2-1 strains) and *Neofusicoccum parvum*.

NA: not applicable - the target was detected in one replicate.

For the three dual labelled probe assays, the PCR mixture contained 1× KAPA Probe Fast Universal qPCR

kit (KK4701) with fungal DNA of different concentrations as template, 300 nM of each of the forward and



**Figure 1.** Standard curves for *Botryosphaeria dothidea* (A), *Neo-fusicoccum mediterraneum* (B) and *N. parvum* (C) specific dual labelled probe (TaqMan) assays. Mean cycle threshold (Ct) values of DNA dilution series for each pathogen were plotted against the log of the DNA concentration (three replicates per concentration of fungal DNA tested). Reactions of undetermined values were set at 40.

reverse primers and probe, and  $1 \times \text{Rox}$  dye. The qPCR conditions were as follows: initial denaturation at 95°C for 3 min, followed by 40 cycles each of 3 sec at 95°C and 30 sec at 60°C, without a final extension period. All reactions were performed in a Step One Plus (Applied Biosystems) real time PCR machine cycler.

As a first step, the quality of all genomic DNA samples was assessed using universal primers ITS4 and ITS5. Amplification was performed with KAPA Taq polymerase using the manufacturer's instructions. PCR conditions consisted of initial denaturation at 95°C for 3 min, followed by 35 cycles each of 30 s at 95°C, 30 s at 55°C, and 60 s at 72°C, and a final extension period of 10 min at 72°C. DNA from the following species was used: *B. dothidea*, *N. parvum*, *N. mediterraneum* MAT1-2-1 strain, while *Diplodia seriata*, *P. chlamydospora* and *N. mediterraneum* MAT1-1-1 strain were included as non-target species.

All three TaqMan assays differentiated the target species from the other two species under study, with no cross-reactions. In particular, the B. dothidea dual labelled probe assay reliably detected target DNA up to 2 pg µL<sup>-1</sup>, at mean Ct 31.3. The N. mediterraneum MAT1-2-1 specific assay detected target DNA up to 20 pg  $\mu$ L<sup>-1</sup>, at mean Ct 31.9. Non-target DNA at 8 ng µL<sup>-1</sup> from *Dip*lodia seriata, P. chlamydospora and N. mediterraneum MAT1-1-1 strain was not detected, while N. parvum and B. dothidea DNA occasionally produced amplification signal at Ct above 35. The N. parvum specific assay detected target DNA up to 2 pg  $\mu$ L<sup>-1</sup> at mean Ct 34.5. Non-target DNA at 8 ng  $\mu$ L<sup>-1</sup> from *Diplodia seriata*, *P*. chlamydospora or N. mediterraneum was not detected, while B. dothidea produced a non-specific signal at Ct 36 (Table 2). For all assays, linear relationships were observed between log of DNA concentrations (serial dilution) and Ct values (Figure 1). Target DNA at 200 pg  $\mu$ L<sup>-1</sup> for all three assays was detected at similar Ct, even if spiked with 8 ng  $\mu$ L<sup>-1</sup> DNA from pistachio (average Ct = 24 for *B. dothidea* and 27 for *N. mediterraneum* and *N.* parvum assays). Pistachio DNA at 8 ng µL<sup>-1</sup> was detected over Ct 35 for all assays.

Due to the high number of species in the Botryosphaeriaceae causing similar symptoms, molecular tools that enable rapid and accurate detection at the pathogen species level are considered to be important for disease management. The present study developed a diagnostic tool to differentiate among B. dothidea, N. mediterraneum and N. parvum. For detection of N. mediterraneum, the assay described here allowed detection only of MAT1-2-1 strains. Another diagnostic tool targeting the MAT1-1-1 gene has yet to be designed. The present study focused on detection of these pathogens because they prevail among species isolated from pistachio, grapevine, citrus (Rumbos and Rumbou, 2001; Vakalounakis et al., 2019; Gusella et al., 2021), and other hosts (including olive, pomegranate, and white willow) (personal communications, Dr Tsopelas; Dr E.J. Paplomatas).

The present study developed three TaqMan assays based on sequence variation of the MAT1-2-1 gene (Lopes *et al.*, 2017) that enabled differentiation between

**Table 3.** List of MAT1.2.1 sequences (Accession number, fungal strain information, host, origin and reference) from species of *Botryosphaeriaceae*, that were aligned using Clustal Omega, in order to detect polymorphisms unique for each of *Botryosphaeria dothidea*, *Neofusicoccum mediterraneum* (MAT1-2-1 strains) and *Neofusicoccum parvum* to design species-specific dual labelled probe qPCR assays.

Fungal strain information		Accession Number	Host	Origin	Reference
1	Botryosphaeria dothidea	Genome database*			
2	Botryosphaeria dothidea MAT1-2-1	OQ632937	Pistacia vera	Greece	This study
3	Neofusicoccum mediterraneum MAT1-2-1	OQ596433	Pistacia vera	Greece	This study
4	Neofusicoccum parvum MAT1-2-1	OQ632936	Pistacia vera	Greece	This study
5	Neofusicoccum mediterraneum strain CAA001 MAT1-2-1	KX505884.1	Pistacia vera	Portugal	Lopes et al., 2017
6	Neofusicoccum parvum strain CMW14085 MAT1-2-1	KX766044.1	N.A.	South Africa	Nagel <i>et al.</i> , 2018
7	Neofusicoccum parvum strain CMW9080 MAT1-2-1	KY612508.1	N.A.	South Africa	Nagel <i>et al.</i> , 2018
8	Neofusicoccum parvum strain CMW9081 MAT1-2-1	KX505872.1	Populus nigra	South Africa	Lopes et al., 2017
9	Neofusicoccum parvum strain CBS 110301 MAT1-2-1	KX505873.1	Populus nigra	The Netherland	s Lopes <i>et al</i> ., 2017
10	Neofusicoccum eucalyptorum strain CAA511 MAT1-2-1	KX505881.1	Eucalyptus globulus	Portugal	Lopes et al., 2017
11	Neofusicoccum mangiferae strain CBS 118531 MAT1-2-1	KX505889.1	Mangifera indica	The Netherland	s Lopes <i>et al.</i> , 2017
12	Neofusicoccum luteum strain CMW9076 MAT1-2-1	KY775143.1	N.A.	South Africa	Nagel <i>et al.</i> , 2018
13	Neofusicoccum algeriense strain CBS 137504 MAT1-2-1	KX505876.1	Vitis vinifera	The Netherlands Lopes et al., 2017	
14	Neofusicoccum algeriense strain CAA322 MAT1-2-1	KX505877.1	Eucalyptus globulus	Portugal	Nagel <i>et al.</i> , 2018
15	Neofusicoccum cordaticola strain CMW14124 MAT1-2-1	KX766043.1	N.A.	South Africa	Nagel <i>et al.</i> , 2018
16	Neofusicoccum cordaticola strain CMW13992 MAT1-2-1	KY612506.1	N.A.	South Africa	Nagel et al., 2018
17	Neofusicoccum kwambonambiense strain CAA755 MAT1-2-1	KX505878.1	Eucalyptus globulus	Portugal	Lopes et al., 2017
18	Neofusicoccum kwambonambiense strain CMW14155 MAT1	KX766045.1	N.A.	South Africa	Nagel et al., 2018
19	Neofusicoccum kwambonambiense strain CMW14023 MAT1	KY612507.1	N.A.	South Africa	Nagel <i>et al.</i> , 2018
20	Neofusicoccum ribis strain CMW7054 MAT1-2-1	KX766041.1	N.A.	South Africa	Nagel <i>et al.</i> , 2018
21	Neofusicoccum ribis strain CMW7772 MAT1-2-1	KY612509.1	N.A.	South Africa	Nagel <i>et al.</i> , 2018
22	Neofusicoccum ribis strain CBS 115475 MAT1-2-1	KX505879.1	Ribes sp.	The Netherland	s Lopes <i>et al</i> ., 2017
23	Neofusicoccum umdonicola strain CMW14106 MAT1-2-1	KX766042.1	N.A.	South Africa	Nagel <i>et al.</i> , 2018
24	Neofusicoccum umdonicola strain MAT1-2-1	KY612510.1	N.A.	South Africa	Nagel <i>et al.</i> , 2018

N.A.: not applicable.

*B. dothidea*, *N. mediterraneum* (MAT1-2-1 strains), and *N. parvum*. The developed diagnostic tool is superior to other differentiation methods for *Botryosphaeriaceae*, as it requires no time-consuming steps such as RLFP-PCR (Slippers *et al.*, 2007) or polyacrylamide electrophoresis for SSCP analyses (Ridgway *et al.*, 2011). Furthermore, its practical application will be important, as it is possible to quantify species within infected plant tissues, potentially contributing to studies of pathogen prevalence and species interactions, and epidemiology of the diseases they cause.

### FUNDING

This Research Project is co-financed by Greece and the European Union (European Social Fund) through the Operational Program "Human Resources Development, Education and Lifelong Learning 2014-2020", and Program EDBM103, "Support for researchers with an emphasis on young researchers-cycle B", in the context of the project "Risk evaluation of development of a resistant form of disease via transfer of mating type genes from the *Botryosphaeria dothidea* species to species of the genus *Neofusicoccum*" (MIS 5048471).

### AVAILABILITY OF DATA AND MATERIAL

Sequence data that support part of the findings of this study are available at the GenBank database under the accession numbers listed in Table 3.

- Aiello D., Gusella G., Fiorenza A., Guarnaccia V., Polizzi G., 2020. Identification of *Neofusicoccum parvum* causing canker and twig blight on ficus carica in Italy. *Phytopathologia Mediterranea* 59: 213–218.
- Batista E., Lopes A., Alves A., 2020. *Botryosphaeriaceae* species on forest trees in Portugal: diversity, distribution and pathogenicity. *European Journal of Plant Pathology* 158: 693–720.
- Batista E., Lopes A., Alves A., 2021. What do we know about *Botryosphaeriaceae*? An overview of a worldwide cured dataset. *Forests* 12.
- Bezerra J.D.P., Crous P.W., Aiello D., Gullino M.L., Polizzi G., Guarnaccia V., 2021. Genetic diversity and pathogenicity of *Botryosphaeriaceae* species associated with symptomatic citrus plants in Europe. *Plants* 10: 1–31.
- Chen S., Li G., Liu F., Michailides T.J., 2015. Novel species of *Botryosphaeriaceae* associated with shoot blight of pistachio. *Mycologia* 107: 780–792.
- Chen S.F., Morgan D.P., Michailides T.J., 2014. *Botryosphaeriaceae* and *Diaporthaceae* associated with panicle and shoot blight of pistachio in California, USA. *Fungal Diversity* 67: 157–179.
- Guarnaccia V., Kraus C., Markakis E., Alves A., Armengol J., ... Gramaje D., 2023. Fungal trunk diseases of fruit trees in Europe: pathogens, spread and future directions. *Phytopathologia Mediterranea*, 563–599 pp.
- Guarnaccia V., Martino I., Tabone G., Brondino L., Gullino M.L., 2020. Fungal pathogens associated with stem blight and dieback of blueberry in northern Italy. *Phytopathologia Mediterranea* 59: 229–245.
- Gusella G., Lawrence D.P., Aiello D., Luo Y., Polizzi G., Michailides T.J., 2021. Etiology of Botryosphaeria Panicle and Shoot Blight of Pistachio (*Pistacia vera*) Caused by *Botryosphaeriaceae* in Italy. *Plant Disease* 106: 1192–1202.
- Kumar S., Stecher G., Li M., Knyaz C., Tamura K., 2018. MEGA X: Molecular Evolutionary Genetics Analysis across Computing Platforms. *Molecular Biology and Evolution* 35: 1547–1549.
- Lazzizera C., Frisullo S., Alves A., Phillips A.J.L., 2008. Morphology, phylogeny and pathogenicity of *Botry-osphaeria* and *Neofusicoccum* species associated with drupe rot of olives in southern Italy. *Plant Pathology* 57: 948–956.
- Lopes A., Phillips A.J.L., Alves A., 2017. Mating type genes in the genus *Neofusicoccum*: Mating strategies and usefulness in species delimitation. *Fungal Biology* 121: 394–404.

- Luna I.J., Besoain X., Saa S., Peach-Fine E., Morales F.C., ... Rolshausen P.E., 2022. Identity and pathogenicity of *Botryosphaeriaceae* and *Diaporthaceae* from Juglans regia in Chile. *Phytopathologia Mediterranea* 61: 79–94.
- Madeira F., Park Y. mi, Lee J., Buso N., Gur T., ... Lopez R., 2019. The EMBL-EBI search and sequence analysis tools APIs in 2019. *Nucleic Acids Research* 47: W636–W641.
- Marsberg A., Kemler M., Jami F., Nagel J.H., Postma-Smidt A., ... Slippers B., 2017. *Botryosphaeria dothidea*: a latent pathogen of global importance to woody plant health. *Molecular Plant Pathology* 18: 477–488.
- Nagel J.H., Wingfield M.J., Slippers B., 2018. Evolution of the mating types and mating strategies in prominent genera in the *Botryosphaeriaceae*. *Fungal Genetics and Biology* 114: 24–33.
- Palavouzis S., Triantafyllopoulou A., Tzima A.K., Paplomatas E.J., 2022. Differentiation between *Botryosphaeria dothidea* and *Neofusicoccum* spp. based on a single nucleotide polymorphism in the ITS region. *Hellenic Plant Protection Journal* 15: 49–56.
- Pavlic D., Slippers B., Coutinho T.A., Wingfield M.J., 2009. Multiple gene genealogies and phenotypic data reveal cryptic species of the *Botryosphaeriaceae*: A case study on the *Neofusicoccum parvum/N. ribis* complex. *Molecular Phylogenetics and Evolution* 51: 259–268.
- Ridgway H.J., Amponsah N.T., Brown D.S., Baskarathevan J., Jones E.E., Jaspers M. V, 2011. Detection of Botryosphaeriaceous species in environmental samples using a multi-species primer pair. *Plant Pathol*ogy 60: 1118–1127.
- Rumbos I., Rumbou A., 2001. Fungi associated with esca and young grapevine decline in Greece. *Phytopathologia Mediterranea* 40: 330–335.
- Slippers B., Smit W.A., Crous P.W., Coutinho T.A., Wingfield B.D., Wingfield M.J., 2007. Taxonomy, phylogeny and identification of *Botryosphaeriaceae* associated with pome and stone fruit trees in South Africa and other regions of the world. *Plant Pathology* 56: 128–139.
- Stempien E., Goddard M.L., Wilhelm K., Tarnus C., Bertsch C., Chong J., 2017. Grapevine *Botryosphaeria* dieback fungi have specific aggressiveness factor repertory involved in wood decay and stilbene metabolization. *PLoS ONE* 12: 1–22.
- Vakalounakis D.J., Ntougias S., Kavroulakis N., Protopapadakis E., 2019. Neofusicoccum parvum and Diaporthe foeniculina associated with twig and shoot blight and branch canker of citrus in Greece. Journal of Phytopathology 167: 527–537.